ORTHO
●●ONE

REF 161 600

KNÖCHEL BANDAGE

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
QUALITY. MADE FOR YOU.

ILLUSTRATIONEN
FIGURES

1

2

3

4


5

SPRACHEN LANGUAGES

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EN English Instructions for use 08

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
Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

INTENDED PURPOSE

This medical device is a support for soft ankle tissue compression consisting of elastic material.

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, upper ankle joint distortion, tendomyopathies.

CONTRAINDICATIONS

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the

region of the body treated, skin diseases in the part of the body treated, latex allergy.

APPLICATION RISKS/IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician.

The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs

- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

FITTING INSTRUCTIONS

- 1** Slip your foot into the slip-on aid provided.
- 2** Then, grasp the upper edge of the support and pull it over the slip-on aid.
- 3** Pull it up further until the tops of the ankles are in the rounded pad recesses. Check the positioning of the support and correct it as applicable.

- 4** When ready, pull the slip-on aid out of the support at the front.

REMOVAL

- 5** To remove the support, please grasp the lower edge and pull it down over the foot.

MATERIAL COMPOSITION

Elastodiene/latex (LA), viscose (CV), polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

CLEANING INFORMATION

-  Delicate wash 30 °C
-  Do not bleach
-  Do not dry in a tumble dryer
-  Do not iron
-  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

USEFUL LIFE/LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link:
www.bort.com/md-eu-contact.

DISPOSAL


Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

DECLARATION OF CONFORMITY

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.ortho-one.net/konformitaet

Status: 10.2024

MD Medical device |  Single patient – multiple use

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BORT GmbH